

ETHICON

MERSILENE*

Braided and Monofilament Polyester Suture



MERSILENE*

Braided and Monofilament Polyester Suture

BRAIDED POLYESTER

Braided

For easy handling and secure knot tying.

Non-Absorbable

Provides permanent tensile strength retention in tissue.

Synthetic

Elicits only minimal tissue reaction. Smooth uniform strand provides easy tissue passage.

Colour

White.

Range

Gauge sizes 5-0 to 2 (U.S.P.)
Supplied as needled sutures, tapes and ligatures.

Sterilisation

By Irradiation.

These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

MONOFILAMENT POLYESTER

Non-Absorbable

Provides prolonged tensile strength retention in tissue.

Monofilament

Smooth surface will not support bacterial growth.

Synthetic

Elicits only minimal tissue reaction.

Knot Security

Higher co-efficient of surface friction than other synthetic monofilament materials gives much improved security.

Tensile Strength

Stronger than other synthetic monofilament suture materials of similar size.

Colour

Green.

Range

Gauge sizes 11-0 and 10-0 (U.S.P.)
Supplied as needled sutures.


Sterilisation


By Irradiation.


These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

MERSILENE*

Braided Polyester Suture



Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting			3-0 2	2-0 3	0 3.5	1 4	2 5	4 6	5 7
 FS-1 24mm	45cm	white	R647H						


Needle length and description	Length	Colour	Gauge						
▼ 1/4 circle SPATULATED			5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
 S-14 8mm	45cm	white	W843	W832					

▼ 1/2 circle SPATULATED			5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5
 S-22 8mm	45cm	white	W833						



MERSILENE*

Monofilament Polyester Suture

Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle Advanced MICRO-POINT spatula			11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG140-6 6.5mm	30cm	green		W1707					
 TG140-6 6.5mm	30cm	green	W1785	W1775					

▼ 1/2 circle Advanced MICRO-POINT spatula			11-0 0.1	10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1
 TG160-6 5.5mm	30cm	green		W1755					

MERSILENE* Tape

Needle length and description	Length & Colour	Width of tape	
<p>48mm</p>			
<p>1/2 circle round bodied heavy double ●</p> <p style="text-align: center;">CTX</p> 	<p>40cm white</p>	<p>5mm</p>	<p>RS22</p>
<p>65mm</p>			
<p>BLUNT POINT 1/2 circle round bodied heavy double ●</p> <p style="text-align: center;">BP-1</p> 	<p>30cm white</p>	<p>5mm</p>	<p>RS21</p>