


**ETHICON**

**PROLENE\***

**Monofilament Polypropylene Suture**



# PROLENE\*

## Monofilament Polypropylene Suture

### Non Absorbable

Provides permanent tensile strength retention in tissue, even in the presence of infection.

### Exceptionally smooth

For easy passage through tissue.

### Monofilament

Smooth surface will not support bacterial growth.

### Synthetic

Elicits only minimal tissue reaction.

### Pliable

For easy knot tying.

### Colour

Blue or clear.

### Range

Gauge sizes 10-0 to 1 (U.S.P.). Supplied as needed sutures.

### Sterilisation

Ethylene Oxide.



As with other synthetic sutures, knot security requires the standard technique of square ties with additional throws as indicated by surgical circumstances and the experience of the operator.














### Caution:

Care should be taken to avoid damaging the surface of the material with surgical instruments, as this could lead to fracture of the material.














These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

# PROLENE\* Monofilament Polypropylene Suture












Needle length and description	Length	Colour	Gauge						
<b>⊙ 1/2 circle ETHIGUARD* Blunt Point</b>			2-0 3	0 3.5	1 4	2 5	4 6	5 7	7 9
 CTB 40mm	100cm	blue			<b>W8998</b>				
 BT-3 45mm	100cm	blue			<b>W8999</b>				

⊙ 3/8 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 BV130-5 6.5mm	13cm	blue	<b>W2775</b>						
 BV130-5 MultiPass* Needles 6.5mm	45cm	blue	<b>W2777</b>						
 BV130-5 6.5mm	60cm	blue	<b>F1890</b>						
 BV175-6 VISI-BLACK* 8mm	60cm	blue	<b>8841H</b>	<b>8335H</b>					
 BV175-6 MultiPass* Needles 8mm	60cm	blue	<b>8741H</b>	<b>8735H</b>					
 BV175-7 VISI-BLACK* 8mm	60cm	blue	<b>8967H</b>	<b>8966H</b>					
 BV175-7 MultiPass* Needles 8mm	60cm	blue	<b>8767H</b>	<b>8766H</b>					
 BV175-8 MultiPass* Needles 9.3mm	60cm	blue		<b>W8702</b>					
 BV175-8 STRAIGHT PACK 9.3mm	4x60cm	blue		<b>W8702S</b>					
 BV-1 VISI-BLACK* 9.3mm	60cm	blue		<b>W8304</b>	<b>W8305</b>				
 BV-1 9.3mm	45cm 60cm	blue blue			<b>W8711</b> <b>W8712</b>				
 BV 11mm	60cm	blue			<b>W8597</b>				
 C-1 VISI-BLACK* 13mm	60cm 75cm	blue blue			<b>W8307</b> <b>W8316</b>	<b>W8321</b>			

# PROLENE\* Monofilament Polypropylene Suture



Needle length and description	Length	Colour	Gauge						
☉ 3/8 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 C-1 13mm	60cm 75cm	blue blue		W8725	W8718 W8706				
 C-1 13mm	75cm	blue				W8890			
☉ 1/2 circle Taper Point			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
 TF-6 6.5mm	75cm	blue		EH7813E					
 RB-2 VISI-BLACK* 13mm	75cm	blue				W8310			
 RB-2 13mm	75cm	blue				W8710			
 RB-1 VISI-BLACK* 17mm	90cm	blue				W8330	W8329		
 RB-1 17mm	90cm	blue				W8556	W8557	W8558	
 SH-2 VISI-BLACK* 20mm	90cm	blue					W8340		
 SH-2 20mm	90cm	blue					W8761		
 SH-1 22mm	75cm	blue						W8770	
 SH VISI-BLACK* 26mm	90cm	blue					W8355	W8354	
 SH 26mm	75cm 90cm 120cm	blue blue blue				F1821	W8521	W8522 W8534	
 JMH-1 VISI-BLACK* 31mm	120cm	blue						X1295H	X1296H










# PROLENE\* Monofilament Polypropylene Suture


Needle length and description		Length	Colour	Gauge					
<b>⊙ 1/2 circle Taper Point</b>				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	3-0 2	2-0 3
MH-1 	31mm	75cm	blue						<b>W295</b>
MH-1 	31mm	90cm	blue					<b>W8525</b>	<b>W8526</b>
<b>⊙ 1/2 circle Taper Point (heavier wire diameters)</b>				2-0 3	0 3.5	1 4	2 5	5 7	7 9
MO-5 	31mm	100cm	blue		<b>W8430</b>				
CT 	40mm	100cm	blue			<b>W742</b>			
MO-2 	40mm	100cm	blue		<b>W8434</b>				
TP-1 	65mm	100cm	blue			<b>W8401</b>			
<b>⊙ Straight Taper Point</b>				2-0 3	0 3.5	1 4	2 5	5 7	7 9
ST-70 	70mm	75cm	blue	<b>W8400</b>					
<b>⊙ 3/8 circle CC Needle</b>				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	3-0 2	2-0 3
CC175-6 VISI-BLACK* 	8mm	60cm	blue	<b>W8101</b>					
CC MultiPass* Needles 	9.3mm	60cm 75cm	blue blue	<b>W8703</b>	<b>W8704</b>	<b>W8807</b> <b>F1832</b>			
CC HEMO-SEAL* MultiPass* Needles 	9.3mm	60cm	blue				<b>W8662</b>		
CC-11 MultiPass* Needles 	11mm	60cm	blue		<b>W8801</b>	<b>W8802</b>			

# PROLENE\*

## Monofilament Polypropylene Suture









Needle length and description		Length	Colour	Gauge						
3/8 circle CC Needle				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	13mm	60cm 75cm	blue blue		W8813	W8707 W8814				
	13mm	75cm	blue				W8721			

1/2 circle CC Needle				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	11mm	60cm	blue				W8803			
	13mm	60cm 75cm	blue blue			W8815	W8816			
	13mm	60cm 75cm	blue blue				W8664	W8665		
	16mm	90cm	blue				W8830	W8831		
	16mm	90cm	blue						W8667	
	20mm	90cm	blue					W8840		
	26mm	90cm 120cm	blue blue					W8845	W8844 W8895	W8843 W8896
	31mm	90cm	blue						W8849	W8850
	40mm	90cm	blue						W8851	W8852



MULTI-CURVE CC Needle				8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
	9.3mm	75cm	blue		EH8021E	EH8025H				


# PROLENE\*

## Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge						
<b>Ⓜ MULTI-CURVE CC Needle</b>			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
ACC-1 VISI-BLACK* 	13mm	75cm	blue			<b>EH8030H</b>	<b>EH8036H</b>		
<b>Ⓜ 3/8 circle TAPERCUT*</b>			8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
V-4 	17mm	75cm	blue			<b>F2829</b>			
V-26 	26mm	75cm 120cm	blue blue				<b>F2845</b> <b>F2850</b>		
<b>Ⓜ 1/2 circle TAPERCUT*</b>			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
V-5 	17mm	90cm	blue	<b>W8935</b>	<b>W8936</b>	<b>W8937</b>			
V-7 	26mm	90cm	blue		<b>W8976</b>	<b>W8977</b>			
V-34 	36mm	100cm	blue				<b>W486</b>		
V-39 	45mm	100cm	blue					<b>W8450</b>	
<b>Ⓜ J Needle TAPERCUT*</b>			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
V-38 	31mm	100cm	blue				<b>W982</b>		

# PROLENE\* Monofilament Polypropylene Suture






Needle length and description	Length	Colour	Gauge						
▼ 1/2 circle Reverse Cutting			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
CP 	40mm	100cm	blue				<b>W8470</b>		
CPX 	48mm	100cm	blue				<b>W8475</b>		

Needle length and description	Length	Colour	Gauge						
▲ 1/2 circle TROCAR POINT			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4	2 5	4 6
TR-30 	31mm	100cm	blue				<b>W8440</b>		










# PROLENE\*



## Monofilament Polypropylene Suture




Needle length and description	Length	Colour	Gauge						
▼ 3/8 circle Reverse Cutting MultiPass* PRIME* Needle			7-0 0.5	6-0 0.7	5-0 1.0	4-0 1.5	3-0 2.0	2-0 3.0	
P-1 	11mm	45cm 75cm	blue blue	<b>MPP2854H</b>	<b>MPP8697H</b>				
P-3 	13mm	45cm 45cm 75cm 75cm	blue clear blue clear		<b>MPP8695H</b>	<b>MPP8698H</b> <b>MPP8605H</b> <b>MPP2831H</b> <b>MPP2858H</b>	<b>MPP8699H</b> <b>MPP8604H</b> <b>MPP2847H</b> <b>MPP2860H</b>		
PS-3 	16mm	45cm 75cm 75cm	blue blue clear		<b>MPP2832H</b>	<b>MPP8681H</b> <b>MPP2833H</b> <b>MPP2859H</b>	<b>MPP7764H</b>		
PS-2 	19mm	45cm 75cm 75cm	blue blue clear			<b>MPP8686H</b>	<b>MPP8682H</b> <b>MPP2834H</b> <b>MPP2861H</b>	<b>MPP8687H</b> <b>MPP2835H</b> <b>MPP2862H</b>	
PS-1 	24mm	45cm	blue				<b>MPP7758H</b>	<b>MPP7757H</b>	<b>MPP7756H</b>

# PROLENE\*





## Monofilament Polypropylene Suture



Needle length and description		Length	Colour	Gauge						
<b>▼ 3/8 circle Reverse Cutting PRIME* Needle</b>				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 P-6	8mm	45cm	blue	<b>W8870T</b>						
 P-1	11mm	45cm	blue	<b>W8868T</b>						
 P-3	13mm	45cm 45cm 75cm	blue clear clear	<b>W8871T</b>	<b>W8872T</b> <b>F2858</b>	<b>8604G</b> <b>F2860</b>				
 PS-3	16mm	45cm 75cm	blue clear		<b>W8882T</b> <b>F2859</b>	<b>W8884T</b>				
 PS-2	19mm	45cm 75cm	blue clear		<b>W8014T</b>	<b>W8015T</b> <b>F2861</b>	<b>F2862</b>			
 PS	26mm	45cm	blue				<b>W8025T</b>	<b>W8026T</b>		
 2 locking buttons	26mm	45cm	blue				<b>W8625</b>			


<b>▼ 3/8 circle Reverse Cutting</b>				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 FS-2	19mm	45cm	blue			<b>W8683</b>				
 FS	26mm	45cm	blue				<b>W8684</b>			

<b>▲ 3/8 circle Conventional Cutting PRIME* Needle</b>				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
 PC-3	16mm	45cm	blue	<b>W8005T</b>	<b>W8006T</b>	<b>W8007T</b>				
 PC-5	19mm	45cm	blue		<b>W8010T</b>	<b>W8011T</b>	<b>W8878T</b>			
 PC-25	26mm	45cm 45cm	blue clear	<b>W8886T</b>	<b>W8881T</b>	<b>W8020T</b> <b>W534</b>	<b>W8021T</b>			










# PROLENE\* Monofilament Polypropylene Suture


Needle length and description		Length	Colour	Gauge					
▲ 3/8 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
FS-3 CONV 	16mm	75cm	blue				<b>W8549</b>		
PCLX  2 beads & 3 collars	36mm	75cm	blue				<b>W626</b>		
PCLX  2 locking buttons	36mm	75cm	blue				<b>W8626</b>		
CEC 135-40 	40mm	100cm	blue					<b>W8689</b>	



▲ 3/8 circle Conventional Cutting (Hand Held)				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
KP-3 CONV  with 40mm retention tubing	90mm	2x50cm	blue						<b>W2989</b>
KP-3 CONV 	90mm	90cm	blue						<b>W987</b>


▲ 1/2 circle Conventional Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	1 4
X-1 CONV 	22mm	45cm	blue				<b>W8731</b>		



# PROLENE\* Monofilament Polypropylene Suture


Needle length and description	Length	Colour	Gauge							
			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5		
<b>▼ Straight Cutting</b>										
 KS-55	55mm	90cm	blue						<b>W8571</b>	
 KS	60mm	75cm	blue						<b>W631</b>	
2 beads & 2 collars  KS	60mm	75cm	blue						<b>W8631</b>	
2 locking buttons  KS	60mm	75cm	blue						<b>W8623</b>	
 KS fine	60mm	75cm	blue					<b>W630</b>		
2 beads & 2 collars  KS fine	60mm	75cm	blue					<b>W8630</b>		
2 locking buttons  KS fine	60mm	75cm	blue					<b>W8621</b>		
 KS fine	60mm	75cm	blue					<b>W621</b>		
 KS-75	75mm	100cm	blue					<b>W751</b>	<b>W752</b>	<b>W753</b>

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
<b>▼ 3/8 circle Advanced MICRO-POINT spatula</b>										
 TG140-6	6.5mm	30cm	blue	<b>W1777</b>						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
<b>▼ 3/8 circle MICRO-POINT spatula</b>										
 TG140-8	6.5mm	23cm	blue	<b>W1711</b>	<b>W1709</b>					
 TG140-8	6.5mm	30cm	blue	<b>W1710</b>						




Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
<b>▼ 1/4 circle TAPERCUT*</b>										
 CIF-4	13mm	20cm	blue	<b>788G</b>						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
<b>▼ Straight MICRO-POINT</b>										
 STC-6	16mm	15cm loop	blue	<b>EH7900G</b>						
 STC-6	16mm	23cm	blue	<b>W1713</b>						

Needle length and description	Length	Colour	Gauge							
			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
<b>▼ 1/4 circle Spatula</b>										
 CTC-6L	13mm	10cm loop	blue	<b>EH7896G</b>						

# PROLENE\*

## Monofilament Polypropylene Suture

Needle length and description	Length	Colour	Gauge							
⊙ 3/8 circle Taper Point			10-0 0.2	9-0 0.3	8-0 0.4	7-0 0.5	6-0 0.7	5-0 1	4-0 1.5	
BV75-3 	3.8mm	13cm	blue	<b>W2790</b>						
BV100-4 	5mm	13cm	blue	<b>W2794</b>	<b>W2783</b>					
BV130-4 	5mm	13cm	blue		<b>W2780</b>					

# PROLENE\*

## Monofilament Polypropylene Suture

### ENDOLOOP\* Ligature

Pre-tied ligatures with Cannula for use during endoscopic surgery

Gauge		blue 45cm
2-0	3	EH464G