



## TROCAR PORTFOLIO

ENDOPATH® XCEL™ Trocars deliver the performance you expect from the market leader

**Abdominal Wall Retention:** Integrated thread design on cannulas of Bladeless, Universal Sleeve, and Dilating Tip Trocars. In Bladeless Trocars, this provides greater stability<sup>1</sup> for minimal trocar slip-outs. Blunt Tip Trocars' adjustable plug with suture stays secures fixation to abdominal wall.

**Seal Durability:** Maintains abdominal insufflation while accommodating a wide range of instrument sizes (4.7mm to 12.9mm).

**Low System Drag Force:** Housing and seal are designed to eliminate instrument hang-ups and provide low peak instrument insertion/extraction force.<sup>1</sup>

**Bladeless Optical Tip:** Direct visualization eliminates blind entry by enabling visualization of tissue layers during insertion. Available in bladeless codes only.



ENDOPATH® XCEL™ Trocars  
Continuing to Put Patients First



Ethicon  
Endo-Surgery

PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

# TROCAR PORTFOLIO FACT SHEET

ENDOPATH® XCEL™ BLADELESS TROCARS		
CODE	DESCRIPTION	LENGTH
B5ST	Stability sleeve, 5mm diameter	75mm
B5SP	Smooth sleeve, 5mm diameter	75mm
B5LT	Stability sleeve, 5mm diameter	100mm
B5LP	Smooth sleeve, 5mm diameter	100mm
B5LPH	Smooth sleeve, 5mm diameter	100mm, handled
B5XT	Stability sleeve, 5mm diameter	150mm
B8LT	Stability sleeve, 8mm diameter	100mm
B11LT	Stability sleeve, 11mm diameter	100mm
B11LP	Smooth sleeve, 11mm diameter	100mm
B11LTH	Stability sleeve, 11mm diameter	100mm, handled
B11LPH	Smooth sleeve, 11mm diameter	100mm, handled
B12SRT	Stability sleeve, 12mm diameter	75mm
B12LT	Stability sleeve, 12mm diameter	100mm
B12LP	Smooth sleeve, 12mm diameter	100mm
B12LTH	Stability sleeve, 12mm diameter	100mm, handled
B12LPH	Smooth sleeve, 12mm diameter	100mm, handled
B12XT	Stability sleeve, 12mm diameter	150mm
B15LT	Stability sleeve, 15mm diameter	100mm

  

ENDOPATH® XCEL™ DILATING TIP TROCARS		
CODE	DESCRIPTION	LENGTH
D5ST	Stability sleeve, 5mm diameter	75mm
D5LT	Stability sleeve, 5mm diameter	100mm
D11LT	Stability sleeve, 11mm diameter	100mm
D12LT	Stability sleeve, 12mm diameter	100mm
D12XT	Stability sleeve, 12mm diameter	150mm

  

ENDOPATH® XCEL™ BLUNT TIP TROCAR		
CODE	DESCRIPTION	LENGTH
H12LP	Smooth sleeve, 12mm diameter	100mm

  

ENDOPATH® XCEL™ UNIVERSAL SLEEVES		
CODE	DESCRIPTION	LENGTH
CB5ST	Stability sleeve, 5mm diameter	75mm
CB5LT	Stability sleeve, 5mm diameter	100mm
CB11LT	Stability sleeve, 11mm diameter	100mm
CB12LT	Stability sleeve, 12mm diameter	100mm

Instruments per Sales Unit: 6

Reference: 1. Compared to ENDOPATH® 5mm and 12mm trocars and/or other competitive trocars. Study EP3-0385 and EP3-0386. Data on file.



**Ethicon**  
**Endo-Surgery**

PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

For answers to product-related questions, call 1-800-USE-ENDO for 24/7 customer support. Or contact your EES Sales Professional.

[www.ethiconendosurgery.com/trocars](http://www.ethiconendosurgery.com/trocars)

© Ethicon Endo-Surgery, Inc. 2010 DSL# 10-1281 ENDO1147